


EMP-7 VET PORTABLE COLOR DOPPLER

VETERINARY DIAGNOSTIC ULTRASOUND SYSTEM



EMP-7 VET PORTABLE COLOR DOPPLER VETERINARY DIAGNOSTIC ULTRASOUND SYSTEM

In nowadays veterinary practice, efficiency and productivity are the most valued for practitioners. EMP-7 Vet can cover all examination needs in your practice; any species and every applications, from small to large animals. With advanced cardiology package to standard everyday scanning, The EMP-7 Vet will be your best choice.



CANINE



FELINE



OVINE



BOVINE



EQUINE




SWINE

Standard Configuration	EMP-7 Vet Main Unit+65C15JK High Frequency Micro Convex Probe
Display	15.6" High Resolution LED Display
Connector	1 probe connectors, 3 probes connector optional
Imaging Modes & Technology	B/2B/4B,M Mode, B+M, Color Doppler, Power Doppler, PW, HPRF,CW
	Tissue Harmonic Imaging
	Tissue Specific Imaging
	Unique Speckle Reduction
	Space Compound Imaging
	Frequency Compound Imaging
	Auto Optimization
	Needle Enhancement for Ultrasound
Optional Configuration	Trapezoid Imaging
	Transducers, Trolley Cart,Gel Warmer

CLINICAL IMAGES




M Mode Bicuspid Valve




Left Ventricular Measurement




Left Ventricle



Portal CFM



Kidney CFM



Carotid Artery

Advanced Measurement & Calculation Packages , Specialized Veterinary Software

Designed Presetting & Body Marks , Advanced Veterinary Abdominal Imaging

Advanced Veterinary Cardiology Imaging


TRANSDUCERS




65C15JK
High Frequency
Micro-convex



65C11JK
Micro-convex



50P10DK
High Frequency
Phased Array



25P20DK
Phased Array



10L25DK
High Frequency
Linear



75L40JK
Linear



35C60VK
Convex



55L60JK
Rectal Linear



EMP

Shenzhen Emperor Electronic Technology Co.,Ltd.

Add:2&3/F, Building 15, No.1008 Songbai Road,

Nanshan District, Shenzhen, 518108, China

Tel:+86-755-36899033 Fax: 86-755-26419886

Email: emp@china-emperor.com

www.china-emperor.com

AUTHORIZED DISTRIBUTOR

