


G30 VET PORTABLE COLOR DOPPLER
VETERINARY DIAGNOSTIC ULTRASOUND SYSTEM



G30 VET COLOR DOPPLER VETERINARY DIAGNOSTIC ULTRASOUND SYSTEM



G30VET is an Color Doppler Veterinary Ultrasound System offering professional software packages, designed presetting and body marks for veterinary use to meet different clinical applications for practiced veterinarians.



CANINE



FELINE



OVINE



BOVINE



SWINE



EQUINE

Standard Configuration	G30 Vet Main Unit+65C15E Micro Convex Probe+Battery+DICOM 3.0.
Display	15" High Resolution LED Display
Connector	1 probe connector
	optional 3 probes connector extension
Weight	6.0kg
Imaging Modes & Technology	B, M, 2B, 4B, B/M, Anatomic M, PDI/DPDI, CFM, PW, CW
	THI Tissue Harmonic Imaging
	TSI Tissue Specific Imaging
	ePure Unique Speckle Reduction
	eSCI Space Compound Imaging
	eFCI Frequency Compound Imaging
	eSpeed Auto Optimization
	eView Panoramic Imaging
Optional Configuration	Trapezoid Imaging
	Mobile Trolley
	Carrying Case

TRANSDUCERS



35C20I
Micro-convex



65C15E
Moco-convex



55L60G
Rectal Linear



10L25K
High Frequency Linear




35C50L
Convex



75L40J
Linear







30P16B
Phased Array



50P10B
Phased Array

CLINICAL IMAGES



Advanced Measurement & Calculation Packages

Specialized Veterinary Software

Designed Presetting & Body Marks



EMP

Shenzhen Emperor Electronic Technology Co.,Ltd.

Add:2&3/F, Building 15, No.1008 Songbai Road,

Nanshan District, Shenzhen, 518108, China

Tel:+86-755-36899033 Fax: 86-755-26419886

Email: emp@china-emperor. com

www.china-emperor.com

AUTHORIZED DISTRIBUTOR

