

G76 VET TROLLEY BASED COLOR DOPPLER
VETERINARY DIAGNOSTIC ULTRASOUND SYSTEM



G76 VET TROLLEY BASED COLOR DOPPLER VETERINARY DIAGNOSTIC ULTRASOUND SYSTEM



In nowadays veterinary practice, efficiency and productivity are most valued for practitioners. G76Vet can cover all examination needs in your practice; any species and every application, from small to large cats and dogs, bovine, equine and to farm animals. With advanced cardiology package to standard everyday scanning, the EMP G76Vet will be your best choice.



CANINE



FELINE



OVINE



BOVINE




EQUINE




SWINE

Standard Configuration	G76 Vet Main Unit+65C15N High frequency Micro-convex+DICOM 3.0.
Display	19"High Resolution LED Display+10.4" Touch Screen
Connector	4 probe connectors
Imaging Modes & Technology	B,M,2B,4B,B/M,PDI/DPDI,CFM,PW,CW, HPRF,3D/4D, Duel Live
	THI Tissue Harmonic Imaging
	TSI Tissue Specific Imaging
	ePure Unique Speckle Reduction
	eSCI Space Compound Imaging
	eFCI Frequency Compound Imaging
	eSpeed Auto Optimization
	eView Panoramic Imaging
Optional Configuration	eHelper Ultrasound Practice Assistant
	Trapezoid Imaging
	Transducers
	ECG, Gel Warmer, Foot Switch

CLINICAL IMAGES




CFM Right Common Carotid Artery



CFM Left kidney



Cardiac Brachial Axis View of Papillary Muscle Section



Cardiac Brachial Pulmonary Artery Reverse Flow



Perfusion of Fetal Heart



B+CFM Left Adrenal Gland

Advanced Measurement & Calculation Packages , Specialized Veterinary Software

Designed Presetting & Body Marks , Advanced Veterinary Abdominal Imaging

Advanced Veterinary Cardiology Imaging

TRANSDUCERS



65C15N
High Frequency
Micro-convex




50P10N
High Frequency
Phased Array



10L25J
High Frequency
Linear



35C20H
Micro Convex



30P16N
Phased Array



75L40N
Linear



35C50N
Convex



35D40J
Volume



55L60H
Rectal Linear



EMP

Shenzhen Emperor Electronic Technology Co.,Ltd.

Add:2&3/F, Building 15, No.1008 Songbai Road,

Nanshan District, Shenzhen, 518108, China

Tel:+86-755-36899033 Fax: 86-755-26419886

Email: emp@china-emperor. com

www.china-emperor.com

AUTHORIZED DISTRIBUTOR

