

*The Innovator of New Ultrasound*

**wisonic**

**Shenzhen Wisonic Medical Technology Co., Ltd.**

Contact address: 6th Floor, NO.6 Building, Pingshan Technology Park,  
Taoyuan Street, Nanshan District, 518055 Shenzhen, P.R.China

Tel : +86 755-8600 7788

Fax: +86 755-8600 7799

Email: [info@wisonic.com](mailto:info@wisonic.com)

Website: [www.wisonic.com](http://www.wisonic.com)

© Shenzhen Wisonic Medical Technology Co., Ltd. All rights reserved. **wisonic** is a registered trademark of Shenzhen Wisonic Medical Technology Co., Ltd. The internal information is for learning only. Shenzhen Wisonic Medical Co., Ltd. reserves the right to modify product specifications and / or terminate the production of any product at any time. Product specifications are subject to change without notice.

H-5900-00-40567 Clivia brochure 8P-V1.0



**Clivia**

Premium POC Ultrasound



**wisonic**

# Versatile and Focused

# Ergonomic Design

With a double-joint articulated arm, Clivia offers optimal viewing angles and comfortable multi-touch gestures



## Bedside▶▶▶

Portable exams in abdominal, cardiac, Ob/Gyn always request excellent mobility and premium image quality. Clivia answers.



## Critical Care▶▶▶

Continuous tracking of each critical patient needs considerable effort to data recording and measuring. Trending graphs and auto-measurement are in demand. Clivia answers.




## Emergency▶▶▶

A fast and confident diagnosis is crucial in every emergency case, which requires quick boot-up, premium image, and excellent portability. Clivia answers.



## General Practice▶▶▶

To satisfy the combined demands of each patient, the device should meet the requirements of multiple departments, including MSK, sports physiotherapy, and cardiac. Clivia answers.



## Clear and Clean

17.3-inch LED HD touch screen, presenting phenomenal imaging quality. With a seamless and waterproof design, Clivia features effortless disinfection.

## Classic and Modern

Clivia delivers dual-workflow, enabling to switch between the classic control panel and fashionable touch-screen arbitrarily.

## Reliable and Convenient

The electric lift provides instant height adjustment to fit your sitting and standing position.

## Versatile and Flexible

Built-in camera supplies remote consultation. 1 to 4 multiprobe adapter expands more applications.

## Safe and Efficient

Fingerprint scanner and password lock.

## Powerful and Durable

Built-in large volume battery and additional battery pack, allow 6+ hours working duration.





# Special and Practical


## Connected and Intelligent

### wiCloud▶▶▶

Clivia equips wiCloud platform, which can implement a variety of cloud applications, such as remote consultation, picture QR code sharing, cloud video tutorial, cloud PACS, cloud quality control, etc.




wiCalc Auto VTI, a valuable tool for in time CO and SV calculation, is very efficient and accurate




Auto Diaphragm Measurement to evaluate diaphragm status and reckon the function of spontaneous respiration



Auto Vascular Measurement, intelligent vascular diameter display, accurate catheter's size recommendation




Holo™ PW Technology, measure blood flow at up to 3 samples simultaneously, assist evaluating vascular stenosis and arteriovenous fistula




wiNerve AI recognition




wiLearn Education Center



wiShow Educational Recording



wiCloud Real-time Remote Consultation



wiNeedle Convex Auto Needle Enhancement



Auto Needle Enhancement, auto needle detection, and auto acoustic beam adjustment


# Full Range of Transducers

A broad range of probes can provide more accurate and comprehensive clinical solutions for our clients. The innovative probe button and probe marking point design, guarantee convenient application within reach.




Transducer Family


## Premium Image Quality




Popliteal Cyst



Carotid Artery Blood Flow




Brachial Plexus with wiNerve




Lower Extremity Blood Flow




Apex 4-chamber



Apex 4-chamber TDI



MV&TV Regurgitation



Renal Blood Flow